

DPD No 1 Soln



The Science of Clean Water

Twenty years of sound scientific research.
The development of proven products.
Offered with wisdom and integrity.
Your guarantee.

Quality is never an accident. It is the
result of high intention, intelligent direction,
committed endeavour and skillful execution.
It is never easy but it is worth it.

46 Bryant Road **P** +64 7 849 7185
PO Box 10333 **F** 0800 725 673
Hamilton, New Zealand **w** sclnz.co.nz



DPD No 1 Soln

At long last the DPD test has been brought into the 21st century. The **DPD No1 Soln** is here. This world first, patent pending breakthrough uses advanced chemistry to combine, concentrate and stabilize all the ingredients.

- You no longer have to wait for a powder or a tablet to dissolve!
- You no longer have the hassle of crushing a tablet!
- You no longer have to carefully measure a powder or two separate liquids!
- All you need do now is add 2 or 3 drops to a comparator or photometer vial or add 1mL to a 100mL titration sample, add your sample and measure! It's instant!

And even more importantly, the test can now be done the way it is supposed to be done; the way it was designed to be done by Dr Palin when he developed the test in the late nineteen fifties. The DPD is put into the vial or flask FIRST and then the sample added.

- Doing it this way, the correct way, produces more accurate, precise, repeatable and reliable results!

And, another plus; it's done much faster.

- The faster the test is done the better as it eliminates interference from combined chlorine compounds and dissolved oxygen.
- Once again producing more accurate, precise, repeatable and reliable results!

And yet another bonus; the NEW, 21st century **DPD No1 Soln** does not go off.

- No need to keep it in a fridge.
- It's stable in ambient conditions for at least two years.
- There is no loss of tablets or powder because of water damage or exposure to air.
- You measure the correct amount every time, time after time.

It is available in two versions; the Eye Dropper and the Titration Soln:

- The 30mL Eye Dropper replaces all the tablets whether for a photometer or a comparator.
- The 200mL Titration Soln replaces the 100g powder and the standard DPD Buffer and DPD Indicator solutions.

Use the ED version for all comparator and photometer testing.

- For a 7mL vial use 2 drops.
- For a 10mL vial use 3 drops.
- The 30mL ED provides 250 tests; the same number of tests as a box of 25 sheets of 10 tablets.

Use the TS version for all titrations.

- Use 1.0mL for a 100mL sample.

COMPARE THE SAVINGS!

COMPARE THE EASE!

YOU'LL SEE A CLEAR ADVANTAGE!

Switch to 21st century technology today. You'll make life so much easier for yourself and you'll save money. From now on the DPD test will be done faster and more accurately than ever before.

And use the **Conc Pot I** and the **Conc Phenol Red** solutions too and save even more time and money!